



ASTM D 1322 - ISO 3014 - IP 57 - DIN 51406

SCOPE:

This test method covers a procedure for determination of the smoke point of kerosene and aviation turbine fuel.

DETERMINATION OF THE SMOKE OF KEROSENE AND AVIATION TURBINE FUEL - REF 941680 -

MAINS CHARACTERISTICS

Smoke point apparatus; complete as per the above mentioned standard methods specifications, delivered with one interchangeable oil container.

DIMENSIONS 200 x 200 x 400 mm - (± 4 kg)

SPARE PARTS AND CONSUMABLES

REF 9416801 Candle (oil container)

REF 9416802 Pack of 12 wicks



NORMALAB FRANCE SAS
14 rue des Lilas - F-76 210 Lintot
Tel. : +332 35 38 59 59 - Fax : +332 35 38 78 55

NORMALAB BELGIUM S.A.
Parc Industriel
28 rue des Pieds-d'Alouette – BE-5100 NANINNE / Belgium
Tel : +32 81 71 19 90 - Fax : +32 81 22 27 01

DISTRIBUTED BY

