

**PROGRAMMABLE OPERATION**

Set time and amplitude for hands free operation

**PULSE MODE**

Prevent heat buildup in temperature sensitive samples

**DIGITAL AMPLITUDE / INTENSITY CONTROL**

Output intensity can be set from 20-100%

**ELAPSED TIME INDICATOR**

Displays duration of sonication

**DISPLAY OF WATTAGE AND JOULES**

Real-time energy monitoring

**OVERLOAD PROTECTION**

Prevents damage to circuitry if a fault occurs

**ROHS COMPLIANT**

Uses lead free components

**PART NO. Q500 INCLUDES:**

- Generator
- Converter
- 1/2" diameter probe
- Power cable
- Converter cable
- Wrench set

Power Rating 500 Watts

Frequency 20 kHz

Programmable Timer 10 Hours

Adjustable Pulse On/Off 1 Second to 1 Minute

Dimensions 8" W x 15.25" L x 8.5" H

Voltage 110V, 50/60Hz

\* 220V not available online, contact us for local distributor