

## Technical Specifications

Specifications	Film Test
Test Range	10~10,000 s/in <sup>2</sup> ·100 mL·1.22 KPa
Pressure Range	0~20 KPa (customization is available for others)
Resolution	0.1 Pa
Specimen Size	Φ12 mm
Test Area	0.025 sq.in. (0.16 cm <sup>2</sup> )
Number of Specimens	3, 2 or 1
Test Gas	O <sub>2</sub> , N <sub>2</sub> , CO <sub>2</sub> and 99.9% dry gas (outside of supply scope)
Gas Supply Pressure	0.4MPa~0.6MPa
Port Size	Φ4 mm PU Tubing
Power Supply	220VAC 50Hz / 120VAC 60Hz
Instrument Dimension	645 mm (L) x 500 mm (W) x 468 mm (H)
Net Weight	60 kg