

**Table 2: Technical Specifications**

<b>Specimen Size</b>	Φ74mm
<b>Specimen Thickness</b>	≤3mm
<b>Test Method</b>	Desiccant Method / Water Method
<b>Standard Test Area</b>	33cm <sup>2</sup>
<b>Carrier Gas</b>	Compressed Air
<b>Drying of Carrier Gas</b>	Long Service-life Desiccator (unnecessary to replace inner core)
<b>Humidifying of Carrier Gas</b>	High-efficiency mist-free humidity generator
<b>Carrier Gas Pressure</b>	≥0.6 MPa
<b>Port Size</b>	Φ6mm PU Tubing