

Technical Specifications

Specifications	MEGA 1500
Load Cell Capacity	500N (Standard) 50N, 100N, 250N, 1000N (Optional)
Accuracy	1% FS
Test Speed	1 ~ 500mm/min (Stepless Speed Change)
Number of Stations	6
Number of Specimens	1 ~ 6
Specimen Width	25mm (Standard Sample Grip)
Clamping Way	Pneumatic Specimen Clamp
Stroke	1200 mm
Gas Supply Pressure	0.6MPa ~ 0.7MPa (Outside of Supply Scope)
Port Size	Φ6mm PU Tubing

Instrument Dimension	850 mm (L) x 700 mm (W) x 2030 mm (H)
Power Supply	220VAC 50Hz / 120VAC 60Hz
Net Weight	168 kg