

Technical Specifications

Specifications	Film Test
Test Range	0.1~100,000 cm ³ /m ² ·24h·0.1MPa (standard)
Number of Specimens	3
Vacuum Resolution	0.1 Pa
Vacuum Degree of Test Chamber	<20 Pa
Test Temperature	15°C~55°C (constant temperature controller is optional)
Accuracy	±0.1°C
Specimen Size	Φ80 mm
Test Area	28.27 cm ²
Test Gas	O ₂ , N ₂ and CO ₂ (outside of supply scope)

Test Pressure	-0.1 MPa ~ +0.1 MPa (standard)
Gas Supply Pressure	0.4 MPa ~ 0.6 MPa
Port Size	Φ6 mm PU Tubing
Instrument Dimension	670 mm (L) x 490 mm (W) x 350 mm (H)
Power Supply	220VAC 50Hz / 120VAC 60Hz
Net Weight	57 kg