

Technical Specifications

Specifications	Film Test
Test Range	0.1~10,000 g/m ² ·24h (standard)
Number of Specimens	1~6 with independent test results
Accuracy	0.01 g/m ² ·24h
Resolution	0.0001 g
Temperature Range	15°C~55°C (standard)
Temperature Accuracy	±0.1°C(standard)
Humidity Range	90%RH~70%RH ^{Note} (standard is 90%RH)
Humidity Accuracy	±1%RH
Air Velocity	0.5~2.5 m/s (customization available)
Specimen Thickness	≤3 mm (customization is available for other thickness)
Test Area	33 cm ²
Specimen Size	Φ74 mm
Test Chamber Size	35 L
Gas Supply	Air
Gas Supply Pressure	0.6 MPa
Port Size	Φ6 mm PU tubing
Instrument Dimension	660 mm (L) x 480 mm (W) x 525 mm (H)
Power Supply	220VAC 50Hz / 120VAC 60Hz
Net Weight	70 kg