

Technical Specifications

Specifications	XLW (EC) Auto Tensile Tester
Load Cell Capacity	500 N (standard) 50 N, 100 N, 250 N (optional)
	750 N, 1000N (Customization Available)
Accuracy	Better than 0.5% FS
Test Speed	Forward 50, 100, 150, 200, 250, 300, 500 mm/min Backward 50, 100, 150, 200, 250, 300, 500 mm/min
Number of Specimens	1
Specimen Width	30 mm (Standard Grip)
	50 mm (Optional Grip)
Clamping Way	Pneumatic Specimen Clamp
Gas Supply	Air(outside of supply scope)
Gas Supply Pressure	0.5 MPa ~ 0.7 MPa
Stroke	950 mm
Instrument Dimension	450 mm (L) x 450 mm (W) x 1410 mm (H)
Power Supply	220VAC 50Hz / 120VAC 60Hz
Net Weight	70 kg