

## Technical Specifications

Specifications	XLW (PC)
Load Cell Capacity	500 N (standard) 50 N, 100 N, 250 N (optional)
	750 N, 1000N (Customization Available)
Accuracy	0.5% FS
Number of Specimens	1
Test Speed	50, 100, 150, 200, 250, 300, 500 mm/min
Specimen Width	30 mm (Standard Grip)
	50 mm (Optional Grip)
Stroke	1000 mm
Instrument Dimension	450 mm (L) x 450 mm (W) x 1410 mm (H)
Power Supply	220VAC 50Hz / 120VAC 60Hz
Net Weight	68 kg