

# Micro-Z CL Specifications

<b>X-ray tube</b>	End window type, Pd anode
<b>Analyzing crystal</b>	RX9 doubly curved crystal
<b>Detector</b>	Sealed proportional counter
<b>Sample chamber</b>	Air
<b>Optical chamber</b>	Vacuum
<b>Power requirements</b>	Single phase 100-120 V 15A 50/60 Hz Single phase 220-240 V 10 A 50/60 Hz
<b>Ambient temperature</b>	15~30°C
<b>Humidity</b>	Less than 75% RH