paralab

novasep
Pure Performance for the Life Sciences

Hipersep®	Pilot	XS	S	M	L	XL	XXL
Productivity range (kg/day) in Batch mode	0.2-4	2-7	2-7	7-33	17-57	29-106	53-165
Productivity range (kg/day) in Cyclojet/Varicol®	8-12	14-21	14-21	66-99	114-171	212-318	330-495
Typical Prochrom® column diameters (mm)	50 80 110 150	150 200	150 200	300 450	450 590	590 800	800 1000
Skid flowrate (L/h)	6*-90	60-200	90-200	150-500	300-1000	500–1600	1000-2500
Gradient accuracy (%)	< 1						
Max pressure (bar)	100						70
Explosion-proof	ATEX zone II Class 1 Divison 2						
Certifications / Equipment standards	CODAP/ASME for columns CE/UL for systems						
Configurable options	High pressure injection module, Eluent & feed filters, Heat exchanger, Backflush configuration valves, pH conductivity, Nitrogen flush for product recovery, Up to 10 collections, Solvent recycling, Cyclojet process						