MiniFlex XpC Specifications

Technique	X-ray diffraction (XRD)
Application benefit	Compact X-ray diffractometer for quality control
Technology	θ-θ X-ray diffractometer equipped with screening/quality control software EasyX
Attributes	• Compact X-ray tube • Compact sample horizontal goniometer • D/teX Ultra250 high-resolution, high-speed 1D X-ray detector
Computer	• Windows 10 Pro 64-bit • Touch panel display (Can be combined with external PC)
Dimensions	700 (W) x 1490 (H) x 890 (D) mm (core unit)
Mass	Approx. 250 kg (core unit)
Power requirements	1φ, 200-230 V 50/60 Hz, 6.5A